



# M40 XS

Ultracompact LED Electric Micromotor



# M40 XS Optic / Non-Optic Electric Micromotor with Internal Water Spray

Inside the NSK M40 XS micromotor is the advanced coreless DC motor technology. DC motor technology delivers smooth, quiet operation, powerful torque and minimal heat generation despite the M40 XS compact size.

## Key Features

- Quiet smooth operation from 60 to 40,000 min<sup>-1</sup>
- Minimal heat generation
- Compact and light coreless electric micromotor
- LED illumination is incorporated inside the optic micromotor (M40 XS)
- Titanium body with Scratch Resistant DURACOAT

## Made Shorter and Lighter for Perfect Balance

M40 micromotor is approximately 10 percent shorter and lighter than previous models, giving it far better weight balance and making it as easy to handle as an air turbine. Its titanium body easily stands up to scratches, wear and tear, and repeated autoclaving.



## Advantages of using - LED light

**High Quality Illumination** : LEDs generate daylight quality illumination which significantly improves visibility when compared to halogen bulbs.

**Durable** : LED light sources last noticeably longer, and are more durable than halogen bulbs.



## Ergonomic Shape

The ergo form design of the M40 XS body offers a seamless connection with the contra-angle attachment handpiece, perhaps even creating the sensation of being the extension of one's hand. This body shape ensures a natural grip, fully supporting even the most detailed surgical procedures.



## Micromotor




Optic MODEL M40 XS ORDER CODE E1135051  
Non-Optic MODEL M40N XS ORDER CODE E1136051

## Connect direct to existing

- Bien-Air® 4VLM tubing
- Bien-Air® 4VR 400 tubing
- Bien-Air® B4VLR tubing

M40 XS →



Bien-Air® is registered trademarks of Bien-Air Dental S.A.



**NAKANISHI INC.**  
www.nsk-inc.com

700 Shimohinata  
Kanuma-shi  
Tochigi 322-8666,  
Japan  
TEL : +81 (0)289-64-3380  
FAX : +81 (0)289-62-5636



**NSK Europe GmbH**  
www.nsk-europe.de

Elly-Beinhorn-Strasse 8  
65760 Eschborn,  
Germany  
TEL : +49 (0)6196 77606-0  
FAX : +49 (0)6196 77606-29

**NSK United Kingdom Ltd**  
www.nsk-uk.com

Office 5, Gateway1000,  
Arlington Business Park,  
Whittle Way,  
Stevenage, SG1 2FP, UK  
TEL : +44 (0)1438 310670  
FAX : +44 (0)1438 310671

**NSK Dental LLC**  
www.nskdental.com

1800 Global Parkway  
Hoffman Estates, IL 60192,  
USA  
TEL : +1 888 675 1675  
FAX : +1 800 838 9328

**NSK Oceania Pty Ltd**  
www.nskoceania.com.au

Unit 22, 198-222 Young St.  
Waterloo, Sydney,  
NSW 2017, Australia  
TEL : +61 (0)2 8306 3000  
FAX : +61 (0)2 8306 3030

**NSK Middle East**  
www.nsk-inc.com

Room 6EA-701, 7th Floor,  
East Wing No.6  
Dubai Airport Free Zone,  
PO Box 54316 Dubai, UAE  
TEL : +971 (0)4-204-5507  
FAX : +971 (0)4-204-5508

**NSK Asia Pte Ltd**  
www.nsk-inc.com

1 Maritime Square,  
#09-33 HarbourFront Centre,  
Singapore 099253  
TEL : +65-6376-5018  
FAX : +65-6376-0233